



During emergencies, the simplest way to evacuate bedridden, unconscious or immobile patients through the staircase is by using the evacuation mat. The evacuation mat will be permanently placed beneath the hospital bed.

## KEY FEATURES

- Lay the mat on the base of the bed.
- Fold the Velcro belt and store it beneath the side of the bed.
- Attach the mat to the mattress using the elastic strap.
- During emergencies, wrap the bedding and secure the patient comfortably and safely with the Velcro strap.
- Strap wings are made up of PVC coated polyester fabric.
- The mat is equipped with a body strap, which is a Velcro belt that secures the body firmly on the bed.

